

---

# Frido Go

## Self Propelled



Shower & Commode Wheelchair



Arcatron Mobility Pvt. Ltd.  
Pune, Maharashtra,  
INDIA

Phone Number  
+91 9049291555

---

## ABOUT US

Our mission is to advance global standards of assisted living and make it dignified and safer. Our motto “Purpose Driven Innovation” steers us to apply engineering and technology for making products aimed at solving unmet needs for focus segments. Products that will live up to our name “Arcatron”, which alludes to ‘accelerating wheels’ and is a not so subtle nod to the Ark of Moses, which possessed unlimited power!

Arcatron is building the next generation of devices to enhance lives of people with limited mobility and elderly care needs. It’s range of shower and commode wheelchairs make users independent, safe and dignified. Also, the products are designed to ease out the tedious, hectic jobs of the care-givers making it safer and hygienic. Now never be worried again about those horrific bathroom visits.

## PROBLEMS FACED BY USERS



### PHYSICAL STRESS

Lifting, shifting and carrying the patient takes a toll on the care-giver




### BACK PROBLEMS

Lifting and carrying causes a lot of back pain. More than 80% care-givers suffer from lower back pain.



### SLIPS AND FALLS

Bathroom slips and falls are among the major causes of injuries among the elderly



### BED RIDDEN

Bed ridden users face the bulk of the problems when it comes to toilet visits



### MENTAL STRESS


It is risky, tough, undignified with the mental stress behind toilet visits and hygiene issues.



### LACK OF SOLUTIONS


Lack of shower and toilet friendly products that are rust and water proof cause difficulty for users

## WHAT OUR PRODUCTS CAN OFFER YOU



### **EASY SHIFTING**

Our products make it easy to shift in and out of the wheelchair painlessly



### **PORTABILITY**

Our products make traveling painless and easy due to the added portability.



### **NO MORE INJURIES**

Our products eliminate injuries that users go through during toilet visits



### **READY IN A MINUTE**

Our products can be assembled and disassembled within 60 seconds without any tools



### **COMMODE FRIENDLY**

Our products are height adjustable making it easy to use commodes with different heights



### **ADJUSTABLE**

Our products adjust to your needs to make it comfortable for prolonged use.



## FRIDO GO SELF PROPELLED

Frido Go is a multi-purpose shower and commode wheelchair which is meant to be the ultimate travel companion. It assembles within 60 seconds without tools, and fits inside a bag. It is designed for the shower with materials built to last and function like none other. With a stainless steel body, swivel full armrests, swivel full footrests, front/side access for hygiene, height adjustable casters, quick release shower big wheels.


- ▶ Made of automotive grade powder coated MS Steel
- ▶ Sturdy construction and stable for people up to 150 KGs.
- ▶ Rolls over the commode directly
- ▶ Removable armrests for easy transfers
- ▶ Medical caster wheels with total lock for maximum safety
- ▶ Waterproof, soft padded PU cushions for maximum seating comfort
- ▶ Tension adjustable backrest for maximum comfort
- ▶ Full footrest which swivels inside for maneuvering in tight spaces
- ▶ Attached commode pan for emergency
- ▶ Removable backrest, armrest and swivel footrest for easy storage
- ▶ Swivel push handles for going further towards the wall or avoiding flush tanks
- ▶ Tension adjustable backrest with a safety belt
- ▶ Mag wheels with quick release and Vajr brakes.




**Armrests** that swivel as well as flip back for an easy shifting in and out



**Brakes** give stability while moving in and out of Frido



**The push handle** can be swiveled in to reach closer to the wall or the flush tank



**Flipping seat cushion** provides comfortable and easy access from all sides




**Adjustable height** for using toilets and shifting over beds of different heights



**Adjustable backrest** tension can be re-tightened as your posture needs dictate



**Folds and packs** into our suitcase within 60 seconds



**The active footrest** ensures smooth turning even in narrow passages



**Built to be resistant** to water and corrosion



**Wheel positions** can be switched so that Frido can be used during toilet visits

ACCESSORIES



**Top seat cover** to use the chair as a regular wheelchair



**Seat belt** for added safety



**Removable bed pan** with cover for bed side use



**Casters** can turn Frido into attendant propelled



**Wrap around bag** makes Frido portable while traveling



**Trolley and big wheel bags** make Frido portable while traveling

## TECHNICAL DETAILS

## FRIDO GO SELF PROPELLED


Narrow model: (N)

Wide model: (W)


Weight capacity: 120kg

Product weight: 19 kg (N) / 19.5 kg (W)

Seat size: 18 x 18 inch



Front



Side



---

[WWW.ARCATRON.COM](http://WWW.ARCATRON.COM)

Fb [@arcatronmobility](#) — Instagram [@myfrido](#) — twitter [@arcatronindia](#) — linkedin [@arcatron-mobility-pvt-ltd](#)